



SALSA Mini TEENS

THE NARROWEST
PERFORMANCE
MID-WHEEL DRIVE
IN THE WORLD!

PERFECT FOR TEENS!



INTRODUCING...

The Zippie Salsa M² Mini Teens, perfect for indoor manoeuvrability and outdoor performance. Perfect for teens the seating system has been specifically designed to accommodate users requiring narrow seat widths with long seat depths. The Salsa M² Mini Teens enable teens to live their lives to the fullest.

Introducing the world's narrowest mid-wheel drive performance powerchair with a chassis width of only 52cm – that's 9cm smaller than our market leading Quickie Salsa M². Unlike other teens indoor powerchairs with the Salsa M² Mini Teens you don't have to compromise on outdoor performance, the Salsa M² Mini Teens uses exactly the same drive base as the Salsa M², just narrower – so all of the great features you get on a bigger outdoor powerchair are packed into the world's narrowest mid-wheel drive powerchair.

INDOOR AGILITY...

The Salsa M² Mini is the narrowest performance mid-wheel drive in the world. This means you can tackle the tightest indoor environments but that's not just it...it's ideal for:

- School buses & wheelchair adapted vehicles – narrow, ultra low seat to floor height and the anti pitch technology not only means getting in and out of WA vehicles is easy, but you'll sit comfortably under the steering wheel and have a stable ride.
- Small Houses & cramped class rooms – with tight floor ways and condensed layouts.
- Through floor lifts – short, narrow and lighter than any other performance mid-wheel drive powerchair – from 100kg means it's ideal for through floor lifts and hoists where other chairs are too big or heavy.

OUTDOOR PERFORMANCE...

Evolved from highly-advanced technology the Salsa M² Mini Teens has proven, patented technology, together with a series of high-performance options.

- Go explore with Zippie's anti-pitch technology eliminates pitching on slopes, and with all your weight located directly above the drive wheel to maximise traction, you've never been in a safer driving position.
- Enjoy a smooth and stable ride, even over rough terrain, with Zippie's patented all-wheel suspension - six independently suspended wheels absorb all those bumps.
- More miles per charge! Take advantage of the ultra-efficient 4mph speed option powered by 36Ah batteries.
- Whatever your performance or practical requirements, the Salsa M² Mini's range of options gives you the optimum balance between an exhilarating drive and great efficiency.



52cm

RIDE IN STYLE...

At Zippie we don't believe in compromise with style that's why the Zippie Salsa M² Teens is our most stylish looking powerchair for teens yet!

The Salsa M² Mini Teens blends Zippie's hallmarks of beautiful sleek design with an modern and sporting stance.

The sleek shroud is available in four different colours (red, blue, white and black) or you can have the exact colour you want through our bespoke Built-4-Me department

The castor arms have been designed with a fresh, modern and slick look and the matt black finish is a great accompaniment for any shroud colour

FUNCTION...

The Salsa M² Mini is so stable that we offer our full range of powered options. That means that you can still choose the options designed for your clinical needs.

- Powered Tilt – if you require pressure redistribution throughout the day then a tilt option can make this easy and comfortable to achieve. It can also be used when hoisting into the chair to guide the pelvis to the back of the seat with ease.
- Powered Lift – is great for raising you up to be at eye-level with your family, friends or colleagues but it can also be used to improve your reach around the house.
- Powered Recline – can often be used to achieve a relax position throughout the day. If you are spending a lot of time in your chair, then a little recline can allow you to change your position every now and again.

TAKE A SEAT...

The Salsa M² Mini is our most dynamic powered wheelchair to date. It made sense that we made it available for all ages.

Teens Seat...

Our teens seat is designed specifically to offer a wide range of adjustments especially to suit teens.

The teens seat offers a range of adjustment which starts at 12" and goes to 16" in width and 16" - 20" in depth. This means that even when you grow, the seat will cater for it and last much longer than the classic teens seating.

Too big or small?

Don't worry the Salsa M² Mini is also available in both kids sizes (12" - 16" width & depth) and adult (16" - 20" width & depth)



SALSA M² MINI // Mid-wheel powered wheelchair

Max. user weight:	75kg (12 stone)	Battery size:	36Ah
Front & rear castor wheels:	7" Solid (18cm)	Max. dynamic stability:	8° (Depending on options)
Drive wheels:	12.5" (32cm)	Overall Chair Width	52cm
Seat width:	12" - 16" (30 - 40cm)	Seat Depth:	16" - 20" (41 - 51cm)
Min turning circle:	43.3" (110cm) - with centre mount	Max. kerb climbing:	3" (7cm)
Max. speed:	2mph and 4mph	Colour choices:	4 colours - red, blue, white & matt black
Max. range (ISO 7176-4):	23km (14 miles)	Transport (ISO 7176-19):	4 point tie-down & Dahl docking system
Battery box to ground clearance:	2.5" (6.5cm)	Chair weight:	97kg - 128kg (Depending on options)

For further information on the full specification, options and accessories, please refer to the order form on our website.

All information is subject to change without notification. Please consult Sunrise Medical with any queries you may have.

If you are visually impaired, this document can be viewed in PDF format at www.SunriseMedical.co.uk or alternatively it is available on request in large text.

